

TRIAD HYDROPHILIC WOUND DRESSING

A prescription is not required

Triad provides a moist wound environment, targets wounds and debrides them, & protects surrounding skin from excess drainage. Ideal for difficult-to-dress wounds with light-to-moderate drainage.

- Remove old dressings and lightly clean wound with Wound Cleanser or soap/water or take a shower
 - Do NOT scrub off hardened cream, just lightly clean
 - You may clean the wound daily (*LOTS of drainage*) to every third day (*NO drainage*). **If outer dressing is wet, CHANGE IT**

- Apply a layer of Triad (**thickness of a dime**) directly on the wound
 - *It is absolutely fine if cream gets on healthy surrounding tissue & if you apply new cream on top of old*

- Cover with an outer dressing, as needed
 - The paste will harden; No outer dressing is needed unless there is drainage.

- The amount of drainage dictates how often the dressing should be changed. We recommend at least every 3rd day, but more often as the dressing is soiled by drainage
 - *Do NOT apply new cream over old cream. Too much build up can create pressure and can cause new wounds*